



# AI in Medical Image Classification

## Warmup

---

8/9 (Sun) 14:00-17:00

中央研究院統計科學研究所  
李易儒博士 森元俊成研究助理

# Warmup

14:00-15:30 學習動機：了解電腦影像與基本Python語法

15:30-17:00 統計模型估計


本課程所使用的肺部X光影像資料來自開源資料庫

<https://www.kaggle.com/paultimothymooney/chest-xray-pneumonia>


# 課前準備

- 下載Anaconda在自己的筆電上

<https://www.anaconda.com/products/individual>



# 在自己的電腦上安裝Python





# 學習動機：了解電腦影像與基本Python語法



- 猜猜看哪一張
- 電腦影像為何？
- 本課程使用的肺部影像介紹(Chest x-ray)
- 如何使用Python 操作影像

猜猜看，哪一張？

# 哪一張是有肺炎的肺部影像？





# 哪一張是有肺炎的肺部影像？





# 哪一張是有肺炎的肺部影像？



# 人工智慧在醫療場域

診斷＝模式(pattern)識別  
有智慧的機器可以幫助醫生…

- 自動分類與檢查影像－增益遠距醫療與急診
- 機器學習輔助第二意見－受益於其他醫師對類似病例的共識
- 自動報告產生－紀錄以影像為本的機器產出
- 診斷細化－簡少多模態訊息
- 自動腫瘤追蹤

# 人如何學習分類？


## 絕對與相對



藍色或綠色？




貓或狗？




蘋果或櫻桃？

# 人如何學習分類？

認知策略 行為選擇




# 人類的分類極限




Tuna  
Chicken  
Goldfish  
Bullfrog  
Catfish  
Tree frog  
Canary  
Cockatiel  
Parakeet  
Elephant  
Owl  
Human


50 Hz-1.1 kHz  
125 Hz-2 kHz  
20 Hz-3 kHz  
100 Hz-3 kHz  
50 Hz-4 kHz  
50 Hz-4 kHz  
250 Hz-8 kHz  
250 Hz-8 kHz  
200 Hz-8.5 kHz  
17 Hz-10.5 kHz  
200 Hz-12 kHz  
31 Hz-19 kHz




The Hearing Range of various animals in terms of Frequency



# 人類的分類極限



# 人類的分類極限？



# 分類：系統性地依據觀察分群入不同類別

## 人

- 吃食物 水
- 經驗累積所歸納的規則與特徵
- 容許無法描述的內隱知識
- 有彈性的邏輯思維

## 電腦


- 吃運算資源
- 可定義的規則
- 知識持續大量累積
- 提供鉅細靡遺的估計資料




# 電腦影像為何？

影像是數字


# 醫學影像作為Data



# 資料數位化使影像可以被電腦分析



# 什麼是醫學影像？



- 包含臨床意義的訊息
- 以輔助治療、診斷、照護為目的



# 醫學影像數位化

醫學資料必需結構化，讓機器可以認得



16 世紀



20 世紀



未來？

# 醫學影像裡的資料



機器端攝影參數



影像本身屬性



影像相關資料

# 量化的資料讓電腦讀得懂 可以分析

## 機器端攝影參數

掃描者姓名  
掃描模組  
掃描序號  
掃描日期  
資料日期  
掃描時間  
掃描時間長  
影像種類  
機器製造商

## 影像本身屬性


掃描序號  
影像厚度  
影像位置  
影像角度  
索引UID  
Slice位置  
欄  
列  
像素Spacing  
Bits Allocated  
Bits Stored  
最小影像像素值  
最大影像像素值  
Window 中央位置  
Window寬  
Rescale截距  
Rescale斜率

## 影像相關資料

姓名  
生日 性別  
年齡 體重  
病歷號  
掃描號  
掃描部位

## 綜合評估與意見

# 醫學影像的格式



- **D**igital **I**maging and **C**ommunications in **M**edicine (DICOM)
- **N**euroimaging **I**nformatics **T**echnology **I**nitiative (NIfTI)
- **N**early **R**aw **R**aster **D**ata (NRRD)



# 醫學影像的格式

## Summary of file formats characteristics

Format	Header	Extension	Data types
Analyze	Fixed-length: 348 byte binary format	.img and .hdr	Unsigned integer (8-bit), signed integer (16-, 32-bit), float (32-, 64-bit), complex (64-bit)
Nifti	Fixed-length: 352 byte binary format <sup>a</sup> (348 byte in the case of data stored as .img and .hdr)	.nii	Signed and unsigned integer (from 8- to 64-bit), float (from 32- to 128-bit), complex (from 64- to 256-bit)
Minc	Extensible binary format	.mnc	Signed and unsigned integer (from 8- to 32-bit), float (32-, 64-bit), complex (32-, 64-bit)
Dicom	Variable length binary format	.dcm	Signed and unsigned integer, (8-, 16-bit; 32-bit only allowed for radiotherapy dose), float not supported

Not all the software support all the specified data types. Dicom, Analyze, and Nifti support color RGB 24-bit; Nifti also supports RGBA 32-bit (RGB plus an alpha-channel)

<sup>a</sup>Nifti has a mechanism to extend the header

# 醫學影像的格式

- Consists of:
  - Metadata  
(eg file size, image dimensions, scan date, patient id)
  - Data (usually big long string of bytes)
- Organized as:
  - Header (metadata) followed by data
  - Containers (each has an identifier, metadata, data)

Lots of information



• DICOM

• HDF5

• MINC


• NIfTI

• Analyze7.5

• RAW



Only data

# 影像是數字




- 最小單位是pixel
- 裡面會有數字
- 黑白影像可以用矩陣表示
- 彩色影像可以用三階張量(RGB)表示

影像是數字



= 16



= 243

# 黑白影像是矩陣



=

0	0	0	0	0	0	0	1	5	8	1	0	0	0
0	0	0	0	0	0	0	1	5	8	1	0	0	0
0	0	0	0	0	120	0	1	5	8	1	0	0	0
0	0	0	0	3	0	0	1	5	8	1	0	0	0
0	0	0	0	0	43	0	1	5	8	1	0	0	0
0	0	0	0	33	0	0	1	5	8	1	0	0	0
0	0	0	0	22	0	0	1	5	8	1	0	0	0
0	0	0	0	0	0	231	165	250	8	1	0	0	0
0	0	0	0	78	0	250	1	254	8	1	0	0	0
0	0	0	112	254	199	222	231	250	231	98	0	0	0
0	33	77	231	250	210	231	250	250	250	86	134	0	0
0	23	89	250	250	228	255	222	222	222	92	169		0
0	56	132	222	222	251	243	255	222	231	132	215	142	0
0	0	255	222	255	222	255	255	222	255	222	0	0	0
0	0	243	255	243	255	255	243	255	243	255	0	0	0

矩陣！


# 本課程使用的肺部影像介紹

## 胸部 X 光影像

# 資料集：胸部 X 光影像 (健康/肺炎)

## 資料來源

<https://www.kaggle.com/paultimothymooney/chest-xray-pneumonia/metadata>




The screenshot displays the Kaggle dataset page for 'Chest X-Ray Images (Pneumonia)' by Paul Mooney. The page features a navigation bar with 'Home', 'Compete', 'Data', 'Notebooks', 'Discuss', 'Courses', and 'More'. The dataset title 'Chest X-Ray Images (Pneumonia)' is prominently displayed, along with the description '5,863 images, 2 categories' and the creator 'Paul Mooney • updated 2 years ago (Version 2)'. A red arrow points to the 'Download (1 GB)' button. The 'Metadata' tab is selected, showing a table with usage information, maintainers, and updates.


Metadata		
Usage Information	License	Other (specified in description)
	Visibility	Public
Maintainers	Dataset owner	Paul Mooney
Updates	Expected update frequency	Not specified
	Last updated	2018-03-25
	Date created	2018-03-22

# 資料集：胸部X光影像 (健康/肺炎)

- 胸部X光影像來自廣州市婦女兒童醫學中心回溯性資料庫
- 1至5歲的小兒科病患影像



健康影像



肺炎影像




# 資料集：胸部 X 光影像 (健康/肺炎)

- 5,863 張X光影像 (JPEG)
- 兩種資料類別(肺炎Pneumonia/正常Normal)
- 資料分別整理在3個資料夾中：train, test, val
- 每個資料夾內分別有兩個類別的影像資料夾

	<b>Normal</b>	<b>Pneumonia</b>
<b>train</b>	<b>1342</b>	<b>3875</b>
<b>test</b>	<b>234</b>	<b>390</b>
<b>val</b>	<b>8</b>	<b>8</b>

# 作業流程




# 聰明的機器，目前有哪些？

- CNN (Wang, et al., 2017; Li et al., 2017)
- CNN + recurrent neural network (Yao et al.,2017)
- DenseNet: as the backbone of CNN (Huang et al.,2017)
- DenseNet-121 with fine tuning : Transfer learning (Rajpurkar, 2017)
- ResNet-50 (Baltruschat et al.,2018)
- VGG 16

# 如何使用 Python 操作影像


利用一張胸部 X 光影像

# 在Anaconda中打開Anaconda Navigator



# 安裝需要的工具

## Environments > base (root) > Open Terminal




# 安裝需要的工具

輸入 conda install **numpy** 後按 enter


輸入 conda install **matplotlib** 後按 enter

輸入 conda install **pillow** 後按 enter

\*部分同學課能會出現已裝好的情形



# 打開Jupyter Notebook





# 下載測試圖片到桌面

[shorturl.at/xLVW2](http://shorturl.at/xLVW2)



test.jpeg

# 開啟新的Jupyter Notebook 點選python 3

The screenshot shows the Jupyter Notebook web interface in a browser window. The address bar indicates the URL is `localhost:8888/tree/Desktop`. The interface includes a top navigation bar with 'Quit' and 'Logout' buttons. Below this, there are tabs for 'Files', 'Running', and 'Clusters'. A message states 'Select items to perform actions on them.' The file browser shows a list of files and folders under the 'Desktop' directory. A 'New' dropdown menu is open, showing options for 'Notebook: Python 3' (with a tooltip 'Create a new notebook with Python 3'), 'Other: Text File', 'Folder', and 'Terminal'. The taskbar at the bottom shows various application icons and the system clock indicating 8:49 PM on 2020/8/7.

Desktop/

localhost:8888/tree/Desktop

Quit Logout

Files Running Clusters

Select items to perform actions on them.

0 / Desktop

Name

Upload New

Notebook: Python 3

Create a new notebook with Python 3

Other: Text File Folder Terminal

3 個月前

9 天前

2 天前

19 小時前

10 個月前

2 個月前

8 個月前

2 個月前

18 天前


3 個月前

1 天前

10 個月前

localhost:8888/tree/Desktop#

下午 08:49 2020/8/7



```
import numpy as np;  
from matplotlib import pyplot;  
from PIL import Image;
```


```
image = Image.open("test.jpeg");  
X = np.array(image);
```

```
print(X);  
print(X.shape);  
pyplot.imshow(X);  
pyplot.imshow(X,cmap='gray');
```

執行：shift+Enter

# 輸入後，畫面呈現如下

jupyter Warmup Last Checkpoint: 12 分鐘前 (unsaved changes)



In [1]:

```
1 import numpy as np;
2 from matplotlib import pyplot;
3 from PIL import Image;
```

科學計算相關的函數庫  
畫圖的函數庫  
讀 jpeg 檔使用

In [2]:

```
1 image = Image.open("test.jpeg");
2 X = np.array(image);
```

讀取剛剛下載的檔案  
將影像轉換成numpy的矩陣

影像是**矩陣**!影像是**矩陣**!影像是**矩陣**!

In [3]:

```
1 print(X);
```

呈現結果的語法

```
[ [ 0 97 92 ... 0 0 1]
  [ 2 98 94 ... 0 0 1]
  [ 2 98 92 ... 0 0 1]
  ...
  [ 0 0 0 ... 0 0 0]
  [ 0 0 0 ... 0 0 0]
  [ 0 0 0 ... 0 0 0]]
```

# 取得影像矩陣Size的語法


```
print (X.shape)
```

```
▶ In [4]: 1 print(X.shape);  
(1484, 1968)
```


# 影像呈現的語法

```
pyplot.imshow(X);  
pyplot.imshow(X,cmap='gray');  
# 以灰階呈現，請輸入cmap='gray'
```

In [5]: 1 pyplot.imshow(X);



In [6]: 1 pyplot.imshow(X,cmap='gray');



# 如何取得特定像素內的數值？

In [7]:

```
1 print(X[0,0]);
```

最左上角  
依矩陣位置編號由零開始

0

In [8]:


```
1 print(X[1,1]);
```

98


0				
	98			



# 如何取得特定像素內的數值？



# 取得特定範圍的像素




In [9]:

```
1 # Show the 左邊100x100  
2 print(X[0:200,0:200])
```

```
[[  0  97  92 ... 152 148 145]  
 [  2  98  94 ... 149 144 140]  
 [  2  98  92 ... 147 142 138]  
 ...  
 [  0  74  64 ... 119 119 119]  
 [  0  74  63 ... 119 119 119]  
 [  0  74  63 ... 118 118 118]]
```


# 取得特定範圍的像素



剪下的圖大小為400\*500  
x軸第 601到第1000個像素  
y軸第1001到第1500個像素


```
In [16]: 1 print(X[600:1000,1000:1500])  
[[194 195 198 ... 87 90 87]  
 [201 200 199 ... 88 90 87]  
 [200 197 194 ... 90 91 87]  
 ...  
 [204 203 203 ... 61 63 65]  
 [203 204 205 ... 61 64 64]  
 [203 204 205 ... 62 65 65]]
```

# 呈現特定範圍的影像



In [10]:

```
1 pyplot.imshow(X[0:800,0:800],cmap='gray');
```



# 建立零矩陣以及所有像素為 1 的矩陣

In [11]:

```
1 # 1484 x 1968 zero matrix
2 A = np.zeros([1484, 1968]);
3 print(A);
```

```
[[0. 0. 0. ... 0. 0. 0.]
 [0. 0. 0. ... 0. 0. 0.]
 [0. 0. 0. ... 0. 0. 0.]
 ...
 [0. 0. 0. ... 0. 0. 0.]
 [0. 0. 0. ... 0. 0. 0.]
 [0. 0. 0. ... 0. 0. 0.]]
```

In [12]:

```
1 B = np.ones([1484, 1968]);
2 print(B);
```

```
[[1. 1. 1. ... 1. 1. 1.]
 [1. 1. 1. ... 1. 1. 1.]
 [1. 1. 1. ... 1. 1. 1.]
 ...
 [1. 1. 1. ... 1. 1. 1.]
 [1. 1. 1. ... 1. 1. 1.]
 [1. 1. 1. ... 1. 1. 1.]]
```

# 影像矩陣與常數的乘法

## 影像矩陣的減法

In [13]:

```
1 C = B * 255;  
2 print(C);
```

```
[[255. 255. 255. ... 255. 255. 255.]  
 [255. 255. 255. ... 255. 255. 255.]  
 [255. 255. 255. ... 255. 255. 255.]  
 ...  
 [255. 255. 255. ... 255. 255. 255.]  
 [255. 255. 255. ... 255. 255. 255.]  
 [255. 255. 255. ... 255. 255. 255.]]
```


In [14]:

```
1 D = C - X;  
2 print(D);
```

```
[[255. 158. 163. ... 255. 255. 254.]  
 [253. 157. 161. ... 255. 255. 254.]  
 [253. 157. 163. ... 255. 255. 254.]  
 ...  
 [255. 255. 255. ... 255. 255. 255.]  
 [255. 255. 255. ... 255. 255. 255.]  
 [255. 255. 255. ... 255. 255. 255.]]
```

## 呈現運算過後的影像矩陣 (反白)

```
In [15]: 1 pyplot.imshow(D,cmap='gray');
```



$$d = 255 - x$$

休息一下